


Info:

tel. 693 465 184

693 465 185

www.dentalholding.com



ioConnect™

Effortless MUA Detection with Your IOS



Open System

Compatible with any dental CAD software for data importation and design.



Effortless MUA Detection with Your IOS

Empower your **existing IOS (Intra-Oral Scanner)** to achieve the same photogrammetry workflow.



Minimized Scanning Area

Focus on scanning **just the blue area** to reduce scanning error. Easy access and shorter scanning time.



All About Options

Multiple sizes (S, M, L) available for you to choose depending on the implant position and the anatomy of each patient.

ioConnect™ Kit

All-on-T (M1.4)

MSRP \$3,000

x 12

x 4

x 4

x 4

x 1

Direct Screw x 20

HS caps

Short/Medium/Long
MUA Detector

1.27 hex
hand driver

Unigrip/Hex1.2/Hex1.27

URIS (M1.6)

MSRP \$3,000

x 12

x 4

x 4

x 4

x 1

Direct Screw x 20

HS caps

Short/Medium/Long
MUA Detector

1.27 hex
hand driver

Unigrip/Hex1.2/Hex1.27

Straumann®

MSRP \$4,000

ST35

x 6

x 2

x 2

x 2

ST46

x 12

x 4

x 4

x 4

x 1

Direct Screw x 20

HS caps

Short/Medium/Long
MUA Detector

Torx
hand driver

Unigrip/Hex1.2/Hex1.27